

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 0.89	
Mean	30.49	Mean	31.81
SD	0.57	SD	0.28
CV	1.86	CV	0.89
Day-Run	Value	Day-Run	Value
2 - 1	29.7	2 - 1	31.5
2 - 2	30.9	2 - 2	31.9
3 - 1	31	3 - 1	32
4 - 1	30	4 - 1	31.5
8 - 1	31.1	8 - 1	32
9 - 1	29.5	9 - 1	31.8
11 - 1	30.9	11 - 1	31.5
12 - 1	29.7	12 - 1	32
15 - 1	30.5	15 - 1	31.9
17 - 1	30.6	17 - 1	31.8
18 - 1	30.9	18 - 1	32.1
21 - 1	29.7	21 - 1	32
22 - 1	31.2	22 - 1	31.5
23 - 1	30.4	23 - 1	32
24 - 1	31	24 - 1	32
25 - 1	30.4	25 - 1	32
28 - 1	31.2	28 - 1	31
29 - 1	30.6	29 - 1	31.9
31 - 1	30	31 - 1	32

Westgard	تفسیر	MCHC	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	5	ماه	1402	سال
Normal Control	Labex	Mean	30.58	SD	0.59	Exp. Date	1402/6/30	Lot. No	02-2305n

