

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	91.83
SD	0.96
CV	1.05
Day-Run	Value
1 - 1	92.2
3 - 1	91.7
4 - 1	92.6
5 - 1	90.6
6 - 1	90.8
7 - 1	91.5
10 - 1	91.5
11 - 1	92.8
12 - 1	92.4
13 - 1	92.4
14 - 1	91.3
15 - 1	90
16 - 1	91.7
17 - 1	93
18 - 1	92
20 - 1	90.5
21 - 1	91.5
22 - 1	93.4
24 - 1	93
25 - 1	92.6
26 - 1	92.5
28 - 1	90
29 - 1	92.2

High - CV : 0.68	
Mean	98.89
SD	0.67
CV	0.68
Day-Run	Value
1 - 1	99.3
3 - 1	99.1
4 - 1	98.9
5 - 1	98.9
6 - 1	98.2
7 - 1	98.5
10 - 1	98.5
11 - 1	98.9
12 - 1	99.1
13 - 1	98.1
14 - 1	98.2
15 - 1	98.6
16 - 1	98.5
17 - 1	100
18 - 1	99.7
21 - 1	97.6
22 - 1	98.9
24 - 1	98
25 - 1	100
26 - 1	100
28 - 1	99
29 - 1	99.5

Westgard	تفسیر	MCV	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	4	ماه	1402	سال
Normal Control		Mean	91.38	SD	1.01	Exp. Date	1402/4/13	Lot. No	

