

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 4.19	
Mean	2.2	Mean	3.9
SD	0.16	SD	0.16
CV	7.42	CV	4.19
Day-Run	Value	Day-Run	Value
1 - 1	2.3	1 - 1	3.9
2 - 1	2.5	2 - 1	3.8
7 - 1	2.2	7 - 1	3.9
8 - 1	2.3	8 - 1	3.7
9 - 1	2.4	9 - 1	3.6
11 - 1	2.2	11 - 1	4
12 - 1	2.2	12 - 1	3.9
14 - 1	2.1	14 - 1	4
15 - 1	2.4	15 - 1	3.9
16 - 1	2.2	16 - 1	4
17 - 1	2.4	17 - 1	3.7
18 - 1	2	18 - 1	4.2
19 - 1	2	19 - 1	3.6
21 - 1	2.3	21 - 1	4.1
22 - 1	2.2	22 - 1	3.9
23 - 1	2	23 - 1	4.1
24 - 1	2	24 - 1	4
25 - 1	2.2	25 - 1	4
26 - 1	2	26 - 1	4.1
28 - 1	1.9	28 - 1	3.8
29 - 1	2.3	29 - 1	3.75
30 - 1	2.3	30 - 1	3.9

Westgard	تفسیر	Magnesium(Mg)	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	2.26	SD	0.14	Exp. Date	1403/1/17	Lot. No	035

