

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 4.4	
Mean	2.26	Mean	3.81
SD	0.14	SD	0.17
CV	6.25	CV	4.4
Day-Run	Value	Day-Run	Value
1 - 1	2.4	1 - 1	3.6
3 - 1	2.2	3 - 1	3.5
4 - 1	2.3	4 - 1	3.9
5 - 1	2.4	5 - 1	3.6
6 - 1	2.4	6 - 1	4
10 - 1	2.5	10 - 1	3.8
11 - 1	2.11	11 - 1	3.7
13 - 1	2.3	13 - 1	3.6
19 - 1	2.1	19 - 1	3.9
20 - 1	2	20 - 1	4.1
21 - 1	2.2	21 - 1	3.8
24 - 1	2.1	24 - 1	4
25 - 1	2.2	25 - 1	3.8
26 - 1	2.3	26 - 1	4
27 - 1	2.3	27 - 1	3.9
28 - 1	2.5	28 - 1	3.8
29 - 1	2.2	29 - 1	3.7
31 - 1	2.2	31 - 1	3.9

Westgard	تفسیر	Magnesium( Mg )	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	2.24	SD	0.13	Exp. Date	1403/1/17	Lot. No	035

