

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	13.6	
2 - 1	13.48	
3 - 1	12.57	
4 - 1	12.75	
6 - 1	13	
7 - 1	13.57	
8 - 1	13.32	
9 - 1	14	
10 - 1	12.95	
11 - 1	13.52	
13 - 1	13.17	
16 - 1	12.25	
17 - 1	12.44	
18 - 1	14.1	
20 - 1	13.75	
21 - 1	13.5	
22 - 1	12.51	
23 - 1	13.5	
24 - 1	13	
25 - 1	12.55	
27 - 1	12.7	
28 - 1	13.4	
29 - 1	13.5	
30 - 1	12.9	
31 - 1	13	

Westgard	تفسیر	Prolactin	آزمایش	SIEMENS Centaur CP	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	13.17	SD	0.52	Exp. Date	1402/10/2	Lot. No	36-02-026
Next Month	Mean	13.16	SD	0.5	CV	3.8			

