

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	12.5	
3 - 1	12.4	
4 - 1	13	
5 - 1	12.5	
6 - 1	12.8	
7 - 1	12.9	
10 - 1	12.9	
11 - 1	12.7	
12 - 1	12.5	
13 - 1	12.9	
14 - 1	12.7	
15 - 1	12.5	
18 - 1	12.53	
19 - 1	12.7	
20 - 1	12.5	
21 - 1	12.5	
24 - 1	12.9	
25 - 1	12.8	
27 - 1	12.7	
28 - 1	12.5	
29 - 1	12.7	
30 - 1	12.5	

Westgard	تفسیر	PT	آزمایش	کواگولومتر آه melung Coagulom	دستگاه	4	ماه	1402	سال
Control	آزمایشگاه	Mean	12.69	SD	0.21	Exp. Date		Lot. No	
Next Month	Mean	12.67	SD	0.18	CV	1.42			

