

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
6 - 1	39.4	
8 - 1	40.33	
9 - 1	39.8	
10 - 1	40.5	
18 - 1	40.1	
22 - 1	39.4	
31 - 1	41.7	
7 - 1	41.55	
9 - 2	40.3	
11 - 1	41.1	
1 - 1	41.2	
1 - 2	41	

Westgard	تفسیر	PTH intact	آزمایش	Euroimmun analyser 1	دستگاه	5	ماه	1402	سال
Control	Monobind	Mean	40.8	SD	0.92	Exp. Date	1402/12/24	Lot. No	
Next Month	Mean	40.53	SD	0.78	CV	1.94			

