

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.86	
Mean	3.45	Mean	8.04
SD	0.11	SD	0.15
CV	3.06	CV	1.86
Day-Run	Value	Day-Run	Value
1 - 1	3.7	1 - 1	8
2 - 1	3.4	2 - 1	8
7 - 1	3.4	7 - 1	8.1
8 - 1	3.5	8 - 1	7.9
9 - 1	3.3	9 - 1	8.4
11 - 1	3.4	11 - 1	8
12 - 1	3.4	12 - 1	7.8
14 - 1	3.3	14 - 1	7.8
15 - 1	3.5	15 - 1	8.07
16 - 1	3.4	16 - 1	7.8
17 - 1	3.4	17 - 1	8.2
18 - 1	3.3	18 - 1	8
19 - 1	3.4	19 - 1	8
21 - 1	3.4	21 - 1	8.2
22 - 1	3.5	22 - 1	8.2
23 - 1	3.5	23 - 1	8
24 - 1	3.5	24 - 1	8
25 - 1	3.6	25 - 1	8.2
26 - 1	3.4	26 - 1	8.2
28 - 1	3.5	28 - 1	8.1
29 - 1	3.6	29 - 1	8
30 - 1	3.6	30 - 1	8

Westgard	تفسیر	Phosphorus (P)	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	3.44	SD	0.12	Exp. Date	1403/1/17	Lot. No	035

