

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.69	
Mean	3.44	Mean	8.01
SD	0.12	SD	0.22
CV	3.48	CV	2.69
Day-Run	Value	Day-Run	Value
1 - 1	3.6	1 - 1	8.2
3 - 1	3.4	3 - 1	8.1
4 - 1	3.3	4 - 1	7.8
5 - 1	3.4	5 - 1	8.4
6 - 1	3.6	6 - 1	8.3
10 - 1	3.4	10 - 1	7.8
11 - 1	3.6	11 - 1	7.9
13 - 1	3.5	13 - 1	8
19 - 1	3.6	19 - 1	8
20 - 1	3.3	20 - 1	8.18
21 - 1	3.6	21 - 1	7.6
24 - 1	3.4	24 - 1	8.1
25 - 1	3.3	25 - 1	7.8
26 - 1	3.4	26 - 1	7.8
27 - 1	3.3	27 - 1	7.8
28 - 1	3.5	28 - 1	8.2
29 - 1	3.3	29 - 1	8.2
31 - 1	3.4	31 - 1	8

Westgard	تفسیر	Phosphorus (P)	آزمایش	BS-800 اتونالایزر	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	3.5	SD	0.15	Exp. Date	1403/1/17	Lot. No	035

