

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

| Normal | | High - CV : 2.23 | |
|---------|-------|------------------|-------|
| Mean | 3.5 | Mean | 8.05 |
| SD | 0.09 | SD | 0.18 |
| CV | 2.7 | CV | 2.23 |
| Day-Run | Value | Day-Run | Value |
| 1 - 1 | 3.4 | 1 - 1 | 8 |
| 2 - 1 | 3.6 | 2 - 1 | 7.8 |
| 3 - 1 | 3.6 | 3 - 1 | 7.7 |
| 6 - 1 | 3.4 | 6 - 1 | 7.9 |
| 8 - 1 | 3.4 | 8 - 1 | 8.1 |
| 9 - 1 | 3.4 | 9 - 1 | 8 |
| 10 - 1 | 3.4 | 10 - 1 | 8 |
| 11 - 1 | 3.5 | 11 - 1 | 7.9 |
| 17 - 1 | 3.4 | 17 - 1 | 8.1 |
| 20 - 1 | 3.5 | 20 - 1 | 8.1 |
| 21 - 1 | 3.5 | 21 - 1 | 8.4 |
| 22 - 1 | 3.5 | 22 - 1 | 8 |
| 23 - 1 | 3.5 | 23 - 1 | 8.2 |
| 24 - 1 | 3.4 | 24 - 1 | 8.2 |
| 25 - 1 | 3.6 | 25 - 1 | 8.3 |
| 27 - 1 | 3.7 | 27 - 1 | 7.9 |
| 28 - 1 | 3.5 | 28 - 1 | 7.9 |
| 29 - 1 | 3.6 | 29 - 1 | 8.2 |
| 30 - 1 | 3.6 | 30 - 1 | 8.3 |
| 31 - 1 | 3.4 | 31 - 1 | 8 |

| | | | | | | | | | |
|----------------|--------|------------------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard | تفسیر | Phosphorus (P) | آزمایش | BS-800 اتوآنالایزر | دستگاه | 3 | ماه | 1402 | سال |
| Normal Control | biorex | Mean | 3.5 | SD | 0.13 | Exp. Date | 1403/1/17 | Lot. No | 035 |

