

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	14.1	
1 - 2	13.23	
3 - 1	12.86	
4 - 1	12.48	
5 - 1	13.2	
6 - 1	12.28	
7 - 1	11.8	
10 - 1	12.61	
11 - 1	13.2	
12 - 1	13.52	
13 - 1	14.17	
14 - 1	12.52	
15 - 1	12.58	
17 - 1	12.89	
18 - 1	14.1	
19 - 1	13.62	
20 - 1	13.08	1:3S
21 - 1	12.56	
22 - 1	12.57	
24 - 1	12.2	
26 - 1	13.7	
27 - 1	12.87	
28 - 1	12.75	
29 - 1	13.52	
31 - 1	12.7	

Westgard	تفسیر	Prolactin	آزمایش	SIEMENS Centaur CP	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	12.98	SD	0.69	Exp. Date	1401/5/1	Lot. No	36-02-028
Next Month	Mean	13	SD	0.62	CV	4.8			

