

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	12.7	
2 - 1	12.8	
3 - 1	13.03	
4 - 1	13.2	
7 - 1	12.7	
8 - 1	14.2	
9 - 1	14	
10 - 1	13.5	
11 - 1	11.95	
14 - 1	12.57	
15 - 1	12	
16 - 1	13.5	
17 - 1	13.7	

Westgard	تفسیر	Prolactin	آزمایش	SIEMENS Centaur CP	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	12.99	SD	0.62	Exp. Date	1402/10/30	Lot. No	36-02-028
Next Month	Mean	13.07	SD	0.7	CV	5.38			

