

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.71	
Mean	7.21	Mean	7.95
SD	0.09	SD	0.14
CV	1.23	CV	1.71
Day-Run	Value	Day-Run	Value
1 - 1	7.2	1 - 1	8
2 - 1	7.2	2 - 1	8
3 - 1	7.3	3 - 1	7.9
6 - 1	7	6 - 1	8.1
8 - 1	7.2	8 - 1	8
9 - 1	7.2	9 - 1	8
10 - 1	7.2	10 - 1	8
11 - 1	7.2	11 - 1	7.7
17 - 1	7	17 - 1	8.2
20 - 1	7.3	20 - 1	7.9
21 - 1	7.3	21 - 1	8
22 - 1	7.2	22 - 1	7.8
23 - 1	7.2	23 - 1	8.1
24 - 1	7.3	24 - 1	7.7
25 - 1	7.2	25 - 1	7.8
27 - 1	7.1	27 - 1	7.8
28 - 1	7.2	28 - 1	8
29 - 1	7.3	29 - 1	7.9
30 - 1	7.2	30 - 1	8
31 - 1	7.3	31 - 1	8.1

Westgard	تفسیر	Protein Total	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Normal Control	biorex	Mean	7.12	SD	0.13	Exp. Date	1403/1/17	Lot. No	035

