

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 2.31	
Mean	7.08	Mean	8.02
SD	0.11	SD	0.19
CV	1.53	CV	2.31
Day-Run	Value	Day-Run	Value
1 - 1	7.2	1 - 1	8
2 - 1	7.2	2 - 1	7.9
7 - 1	7.3	7 - 1	8.2
8 - 1	7.3	8 - 1	7.9
9 - 1	7	9 - 1	7.8
11 - 1	7	11 - 1	7.9
11 - 1	7	12 - 1	7.8
12 - 1	7.2	12 - 2	8.1
14 - 1	7	14 - 1	7.9
15 - 1	7	15 - 1	7.8
16 - 1	7	16 - 1	8
17 - 1	7.1	17 - 1	8
18 - 1	7	18 - 1	7.9
19 - 1	7.1	21 - 1	7.9
21 - 1	7	22 - 1	8
22 - 1	7	23 - 1	8.3
23 - 1	7.1	24 - 1	7.8
24 - 1	7.2	25 - 1	8.4
25 - 1	7.1	26 - 1	8.3
26 - 1	7	26 - 2	8.3
28 - 1	7	28 - 1	7.9
29 - 1	7.1	29 - 1	8.2
30 - 1	6.9	30 - 1	8.1

Westgard	تفسیر	Protein Total	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	7.19	SD	0.15	Exp. Date	1403/1/17	Lot. No	035

