

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.28	
Mean	4.35	Mean	5.5
SD	0.05	SD	0.07
CV	1.2	CV	1.28
Day-Run	Value	Day-Run	Value
2 - 1	4.33	2 - 1	5.45
2 - 2	4.3	2 - 2	5.55
3 - 1	4.39	3 - 1	5.45
4 - 1	4.38	4 - 1	5.55
8 - 1	4.42	8 - 1	5.52
9 - 1	4.3	9 - 1	5.51
11 - 1	4.31	11 - 1	5.45
12 - 1	4.42	12 - 1	5.46
15 - 1	4.31	15 - 1	5.54
17 - 1	4.36	17 - 1	5.52
18 - 1	4.29	18 - 1	5.54
21 - 1	4.4	21 - 1	5.45
22 - 1	4.3	22 - 1	5.37
23 - 1	4.39	23 - 1	5.6
24 - 1	4.25	24 - 1	5.64
25 - 1	4.4	25 - 1	5.38
28 - 1	4.34	28 - 1	5.58
29 - 1	4.33	29 - 1	5.5
31 - 1	4.42	31 - 1	5.46

Westgard	تفسیر	RBC	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	5	ماه	1402	سال
Normal Control	Labex	Mean	4.32	SD	0.06	Exp. Date	1402/6/30	Lot. No	02-2305n

