

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.68	
Mean	0.98	Mean	4.29
SD	0.06	SD	0.12
CV	6.58	CV	2.68
Day-Run	Value	Day-Run	Value
1 - 1	1	1 - 1	4.25
3 - 1	0.92	3 - 1	4.42
4 - 1	1	4 - 1	4.11
5 - 1	0.93	5 - 1	4.41
6 - 1	1.05	6 - 1	4.21
10 - 1	1.03	10 - 1	4.38
11 - 1	0.9	11 - 1	4.38
13 - 1	1.04	13 - 1	4.11
19 - 1	1.05	19 - 1	4.3
20 - 1	0.92	20 - 1	4.25
21 - 1	1.03	21 - 1	4.46
24 - 1	1.1	24 - 1	4.27
25 - 1	0.91	25 - 1	4.11
26 - 1	0.92	26 - 1	4.3
27 - 1	0.92	26 - 2	4.3
28 - 1	1.04	27 - 1	4.47
29 - 1	0.93	28 - 1	4.19
29 - 2	1	29 - 1	4.4
31 - 1	0.9	31 - 1	4.26

Westgard	تفسیر	Bili T	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Normal Control	biorex	Mean	0.97	SD	0.05	Exp. Date	1403/1/17	Lot. No	035

