

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	0.97
SD	0.05
CV	5.34
Day-Run	Value
1 - 1	1
2 - 1	1
3 - 1	0.92
6 - 1	0.9
8 - 1	1
9 - 1	0.92
10 - 1	1.02
11 - 1	0.95
17 - 1	1.03
20 - 1	0.9
21 - 1	1.03
22 - 1	0.9
23 - 1	1
24 - 1	0.98
25 - 1	0.9
27 - 1	1.03
28 - 1	0.93
29 - 1	1.04
29 - 1	1.04
30 - 1	0.95
31 - 1	0.96

High - CV : 3.33	
Mean	4.29
SD	0.14
CV	3.33
Day-Run	Value
1 - 1	4.35
2 - 1	4.4
3 - 1	4.42
6 - 1	4.1
8 - 1	4.3
9 - 1	4.33
10 - 1	4.02
11 - 1	4.3
17 - 1	4.15
20 - 1	4.4
21 - 1	4.11
22 - 1	4.5
23 - 1	4.1
24 - 1	4.44
25 - 1	4.35
27 - 1	4.3
28 - 1	4.05
29 - 1	4.4
29 - 2	4.4
30 - 1	4.28
31 - 1	4.4

Westgard	تفسیر	Bili T	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Normal Control	biorex	Mean	0.95	SD	0.04	Exp. Date	1403/1/17	Lot. No	035

