

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	10.5	
2 - 1	10.52	
3 - 1	10.7	
4 - 1	11.52	
7 - 1	11	
8 - 1	10.8	
9 - 1	12	
10 - 1	11.8	Cusum 3S
11 - 1	10.8	
14 - 1	10	
15 - 1	11.1	
16 - 1	12	
17 - 1	10.7	
18 - 1	11.7	
19 - 1	11	
21 - 1	10.9	
22 - 1	10.7	
23 - 1	11	
24 - 1	12	
25 - 1	11	
26 - 1	10.7	
28 - 1	11.2	
29 - 1	11.5	
30 - 1	11.2	
31 - 1	10.53	

Westgard	تفسیر	T4	آزمایش	SIEMENS Centaur CP	دستگاه	5	ماه	1402	سال
Control	biorex	Mean	10.78	SD	0.64	Exp. Date	1402/10/30	Lot. No	36-02-028
Next Month	Mean	11.07	SD	0.53	CV	4.79			

