

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	299	
3 - 1	311.4	
4 - 1	302.2	
5 - 1	299.5	
6 - 1	316	
10 - 1	301	
11 - 1	315	
13 - 1	302.7	
19 - 1	309.4	
20 - 1	319.9	1:2S
21 - 1	303.6	
24 - 1	311	
25 - 1	319.6	1:2S
26 - 1	314	
27 - 1	310.1	
28 - 1	292.1	
29 - 1	316.1	
31 - 1	325.6	1:2S

Westgard	تفسیر	TIBC	آزمایش	اتوآنالایزر BS-800	دستگاه	4	ماه	1402	سال
Control	biorex	Mean	305.58	SD	6.74	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	309.34	SD	8.85	CV	2.86			

