

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	298.6	
2 - 1	311	
3 - 1	313	
6 - 1	296.2	
8 - 1	308.5	
9 - 1	312	
10 - 1	299	
11 - 1	308.4	
17 - 1	311	
20 - 1	292	
21 - 1	296.3	
22 - 1	305	
23 - 1	311.1	
24 - 1	296.7	
25 - 1	304.7	
27 - 1	311	
28 - 1	308.9	
29 - 1	307.1	
30 - 1	306	
31 - 1	315	1:2S

Westgard	تفسیر	TIBC	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Control	biorex	Mean	302.8	SD	5.72	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	305.58	SD	6.74	CV	2.21			

