

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

| Day-Run | Value | Description       |
|---------|-------|-------------------|
| 1 - 1   | 317.2 |                   |
| 2 - 1   | 323   |                   |
| 7 - 1   | 334.3 | 1:2S              |
| 8 - 1   | 318   |                   |
| 9 - 1   | 314   |                   |
| 11 - 1  | 317   |                   |
| 12 - 1  | 298.1 |                   |
| 14 - 1  | 320   |                   |
| 15 - 1  | 323.5 |                   |
| 16 - 1  | 320.7 |                   |
| 17 - 1  | 307   |                   |
| 18 - 1  | 314   |                   |
| 19 - 1  | 317   |                   |
| 21 - 1  | 296   |                   |
| 22 - 1  | 327   |                   |
| 23 - 1  | 311.2 | 1:2S              |
| 23 - 1  | 311.2 | 1:2S 2:2S         |
| 24 - 1  | 327.7 | 1:2S 4:1S         |
| 25 - 1  | 305   | 1:2S 2:2S<br>4:1S |
| 26 - 1  | 298.6 |                   |
| 28 - 1  | 318.5 |                   |
| 29 - 1  | 311.4 |                   |
| 30 - 1  | 311.2 |                   |

|            |        |        |        |                    |        |           |           |         |     |
|------------|--------|--------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard   | تفسیر  | TIBC   | آزمایش | اتوآنالایزر BS-800 | دستگاه | 5         | ماه       | 1402    | سال |
| Control    | biorex | Mean   | 309.34 | SD                 | 8.85   | Exp. Date | 1403/1/17 | Lot. No | 035 |
| Next Month | Mean   | 315.84 | SD     | 10.39              | CV     | 3.29      |           |         |     |

