

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	173.41
SD	2.04
CV	1.18
Day-Run	Value
1 - 1	172
2 - 1	173
7 - 1	178
8 - 1	174
9 - 1	170
11 - 1	173
12 - 1	175
14 - 1	172
15 - 1	174
16 - 1	173
17 - 1	176
18 - 1	172
19 - 1	175
21 - 1	170
22 - 1	171
23 - 1	173
24 - 1	172
25 - 1	175
26 - 1	173
28 - 1	174
29 - 1	173
30 - 1	177

High - CV : 1.7	
Mean	228.68
SD	3.9
CV	1.7
Day-Run	Value
1 - 1	227
2 - 1	229
7 - 1	221
8 - 1	222
9 - 1	233
11 - 1	232
12 - 1	233
14 - 1	225
15 - 1	232
16 - 1	224
17 - 1	230
18 - 1	230
19 - 1	235
21 - 1	225
22 - 1	231
23 - 1	224
24 - 1	229
25 - 1	233
26 - 1	226
28 - 1	229
29 - 1	231
30 - 1	230

Westgard	تفسیر	Triglycerides	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	172.44	SD	2.33	Exp. Date	1403/1/17	Lot. No	035

