

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

| Normal  |       |
|---------|-------|
| Mean    | 33.99 |
| SD      | 1.26  |
| CV      | 3.69  |
| Day-Run | Value |
| 1 - 1   | 35    |
| 3 - 1   | 35.2  |
| 4 - 1   | 33.7  |
| 5 - 1   | 34.3  |
| 6 - 1   | 33.5  |
| 10 - 1  | 33.7  |
| 11 - 1  | 33.3  |
| 13 - 1  | 33.4  |
| 19 - 1  | 33.3  |
| 20 - 1  | 35    |
| 21 - 1  | 36.6  |
| 24 - 1  | 36    |
| 25 - 1  | 33.4  |
| 26 - 1  | 33    |
| 27 - 1  | 32    |
| 28 - 1  | 35    |
| 29 - 1  | 33.4  |
| 31 - 1  | 32    |

| High - CV : 2.7 |        |
|-----------------|--------|
| Mean            | 134.19 |
| SD              | 3.62   |
| CV              | 2.7    |
| Day-Run         | Value  |
| 1 - 1           | 136.1  |
| 3 - 1           | 137    |
| 4 - 1           | 134.8  |
| 5 - 1           | 140    |
| 6 - 1           | 139    |
| 10 - 1          | 129.5  |
| 11 - 1          | 140    |
| 13 - 1          | 128.9  |
| 19 - 1          | 136.2  |
| 20 - 1          | 133.7  |
| 21 - 1          | 137    |
| 24 - 1          | 130.4  |
| 25 - 1          | 130    |
| 26 - 1          | 134    |
| 27 - 1          | 131.6  |
| 28 - 1          | 135    |
| 29 - 1          | 130.9  |
| 31 - 1          | 131.3  |

|                |        |      |        |                    |        |           |           |         |     |
|----------------|--------|------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard       | تفسیر  | Urea | آزمایش | BS-800 اتوآنالایزر | دستگاه | 4         | ماه       | 1402    | سال |
| Normal Control | biorex | Mean | 34.58  | SD                 | 0.93   | Exp. Date | 1403/1/17 | Lot. No | 035 |

