

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	34.01
SD	1.02
CV	3.01
Day-Run	Value
1 - 1	34.7
2 - 1	33.2
7 - 1	34.2
8 - 1	33.1
9 - 1	32
11 - 1	33.9
12 - 1	35
14 - 1	34.5
15 - 1	33
16 - 1	33.9
17 - 1	36
18 - 1	34.2
19 - 1	33.5
21 - 1	33.4
22 - 1	33.7
23 - 1	34.5
24 - 1	34.2
25 - 1	32.8
26 - 1	33.7
28 - 1	36
29 - 1	33.2
30 - 1	35.6

High - CV : 1.65	
Mean	134.31
SD	2.22
CV	1.65
Day-Run	Value
1 - 1	132.3
2 - 1	133.8
7 - 1	135.7
8 - 1	133.5
9 - 1	133.1
11 - 1	135.8
12 - 1	132
14 - 1	134.2
15 - 1	135
16 - 1	136
17 - 1	140
18 - 1	132
19 - 1	135.8
21 - 1	133.8
22 - 1	133
23 - 1	137.5
24 - 1	135.1
25 - 1	135
26 - 1	131
28 - 1	135
29 - 1	130
30 - 1	135.3

Westgard	تفسیر	Urea	آزمایش	BS-800 اتوآنالایزر	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	33.99	SD	1.26	Exp. Date	1403/1/17	Lot. No	035

