

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 3.08	
Mean	4.46	Mean	9.28
SD	0.11	SD	0.29
CV	2.45	CV	3.08
Day-Run	Value	Day-Run	Value
1 - 1	4.7	1 - 1	9.6
2 - 1	4.4	2 - 1	9.5
3 - 1	4.6	3 - 1	9.6
6 - 1	4.5	6 - 1	9.4
8 - 1	4.5	8 - 1	9.5
9 - 1	4.5	9 - 1	9.4
10 - 1	4.5	10 - 1	9.4
11 - 1	4.5	11 - 1	9.3
20 - 1	4.4	17 - 1	9.6
21 - 1	4.4	20 - 1	8.9
22 - 1	4.4	21 - 1	8.9
23 - 1	4.5	22 - 1	9.7
24 - 1	4.4	23 - 1	9
25 - 1	4.3	24 - 1	9
27 - 1	4.6	25 - 1	8.9
28 - 1	4.3	27 - 1	9.1
30 - 1	4.5	30 - 1	8.9
31 - 1	4.3	31 - 1	9.3

Westgard	تفسیر	Uric Acid	آزمایش	اتوآنالایزر BS-800	دستگاه	3	ماه	1402	سال
Normal Control	biorex	Mean	4.5	SD	0.12	Exp. Date	1403/1/17	Lot. No	035

