

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.63	
Mean	4.26	Mean	8.96
SD	0.07	SD	0.15
CV	1.56	CV	1.63
Day-Run	Value	Day-Run	Value
1 - 1	4.2	1 - 1	8.9
2 - 1	4.3	2 - 1	8.9
7 - 1	4.2	7 - 1	9.1
8 - 1	4.2	8 - 1	9.1
9 - 1	4.2	9 - 1	9
11 - 1	4.3	11 - 1	9.1
12 - 1	4.3	12 - 1	8.6
14 - 1	4.2	14 - 1	9.1
15 - 1	4.3	15 - 1	8.9
16 - 1	4.2	16 - 1	9.1
17 - 1	4.3	17 - 1	9
18 - 1	4.3	18 - 1	9
19 - 1	4.4	19 - 1	9.2
21 - 1	4.2	21 - 1	9
22 - 1	4.3	22 - 1	8.8
23 - 1	4.2	23 - 1	8.9
24 - 1	4.2	24 - 1	8.8
25 - 1	4.2	25 - 1	8.9
26 - 1	4.3	26 - 1	8.8
28 - 1	4.3	28 - 1	8.9
29 - 1	4.4	29 - 1	9.2
30 - 1	4.2	30 - 1	8.9

Westgard	تفسیر	Uric Acid	آزمایش	اتوآنالایزر BS-800	دستگاه	5	ماه	1402	سال
Normal Control	biorex	Mean	4.3	SD	0.07	Exp. Date	1403/1/17	Lot. No	035

