

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	25	
3 - 1	25.2	
4 - 1	24.6	
5 - 1	24.7	
6 - 1	24.3	
10 - 1	24.9	
11 - 1	24.8	
12 - 1	25.1	
13 - 1	24.3	
14 - 1	25	
15 - 1	24.4	
17 - 1	24.5	
18 - 1	24.8	
19 - 1	25	
20 - 1	24.6	
22 - 1	24.8	
24 - 1	25.2	
25 - 1	24.6	
2 - 1	25.2	
2 - 2	24.8	
2 - 3	24.9	
2 - 4	24.5	

Westgard	تفسیر	V D3	آزمایش	HPLC Agilent 1100 -1	دستگاه	4	ماه	1402	سال
Control	آزمایشگاه	Mean	24.75	SD	0.51	Exp. Date	1402/4/31	Lot. No	
Next Month	Mean	24.78	SD	0.28	CV	1.14			

