

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	29.2	
2 - 2	28.1	
3 - 1	28.7	
4 - 1	30.6	
7 - 1	27.4	
8 - 1	29.7	
9 - 1	28.5	
11 - 2	30	
14 - 1	30.8	
15 - 1	29.5	
16 - 1	28.6	
17 - 1	30.5	
18 - 1	31	
19 - 1	30.1	
21 - 1	29.5	
22 - 1	31.4	
23 - 1	29.5	
25 - 1	28.9	
26 - 1	29	
28 - 1	30	
29 - 1	27.7	
30 - 1	30.6	
31 - 1	31	

Westgard	تفسیر	V D3	آزمایش	HPLC Agilent 1100 -1	دستگاه	5	ماه	1402	سال
Control	آزمایشگاه	Mean	29.84	SD	2.47	Exp. Date	1403/6/5	Lot. No	
Next Month	Mean	29.58	SD	1.1	CV	3.73			

