

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	6.73
SD	0.15
CV	2.25
Day-Run	Value
1 - 1	6.7
4 - 1	6.8
5 - 1	6.6
6 - 1	6.8
7 - 1	6.9
10 - 1	6.7
11 - 1	6.7
12 - 1	6.5
13 - 1	6.6
14 - 1	6.6
15 - 1	6.8
17 - 1	6.9
18 - 1	7
19 - 1	7
20 - 1	6.9
20 - 2	6.5
22 - 1	6.5
22 - 2	6.6
25 - 1	6.6
26 - 1	6.7
27 - 1	6.8
28 - 1	6.8
29 - 1	6.7

High - CV : 1.94	
Mean	17.82
SD	0.35
CV	1.94
Day-Run	Value
3 - 1	18.2
4 - 1	18.3
5 - 1	17.9
6 - 1	18
7 - 1	17.8
10 - 1	17.7
11 - 1	17.6
12 - 1	17.5
13 - 1	17.9
14 - 1	18.6
15 - 1	17.6
17 - 1	18.1
18 - 1	18.1
19 - 1	18
20 - 1	17.6
22 - 1	17.1
22 - 2	17.2
25 - 1	17.8
26 - 1	17.5
27 - 1	17.9
28 - 1	17.7
29 - 1	17.9

Westgard	تفسیر	WBC	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	4	ماه	1402	سال
Normal Control		Mean	6.73	SD	0.15	Exp. Date	1402/5/11	Lot. No	

