

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.53	
Mean	6.67	Mean	17.57
SD	0.13	SD	0.27
CV	1.99	CV	1.53
Day-Run	Value	Day-Run	Value
2 - 1	6.7	2 - 1	17.3
2 - 2	6.8	2 - 2	17.5
3 - 1	6.8	3 - 1	17.8
4 - 1	6.8	4 - 1	17.4
8 - 1	6.9	8 - 1	17.2
9 - 1	6.7	9 - 1	17.5
11 - 1	6.6	11 - 1	18
12 - 1	6.5	12 - 1	17.5
15 - 1	6.7	15 - 1	17.5
17 - 1	6.8	17 - 1	17.9
18 - 1	6.5	18 - 1	17.2
21 - 1	6.5	21 - 1	17.5
22 - 1	6.7	22 - 1	17.5
23 - 1	6.85	23 - 1	17.7
24 - 1	6.6	24 - 1	18
25 - 1	6.5	25 - 1	17.7
23 - 2	6.6	23 - 2	17.5
29 - 1	6.5	29 - 1	18
31 - 1	6.7	31 - 1	17.2

Westgard	تفسیر	WBC	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	5	ماه	1402	سال
Normal Control	Labex	Mean	6.73	SD	0.15	Exp. Date	1402/6/30	Lot. No	02-2305n

