

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
3 - 1	69.7	
5 - 1	69.8	
7 - 1	69.77	
11 - 1	70.03	
13 - 1	69	
14 - 1	68.8	
17 - 1	70.02	
18 - 1	71	
20 - 1	70.06	
24 - 1	69.1	
25 - 1	69.88	
26 - 1	68.6	
28 - 1	69.8	
31 - 1	70.1	

Westgard	تفسیر	Selenium	آزمایش	Atomic Absorption Spectrometry	دستگاه	4	ماه	1402	سال
Control	آزمایشگاه	Mean	69.88	SD	0.7	Exp. Date		Lot. No	
Next Month	Mean	69.69	SD	0.63	CV	0.9			

