

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.87	
Mean	28.59	Mean	135.74
SD	0.76	SD	2.54
CV	2.66	CV	1.87
Day-Run	Value	Day-Run	Value
3 - 1	29.5	3 - 1	134.2
4 - 1	28	4 - 1	135.1
5 - 1	29.9	5 - 1	133.7
6 - 1	29	6 - 1	136.8
8 - 1	29.2	8 - 1	132.2
9 - 1	28	9 - 1	138.1
10 - 1	29	10 - 1	133
12 - 1	28	12 - 1	139
13 - 1	28.2	13 - 1	137.7
15 - 1	28.5	15 - 1	131.5
16 - 1	27.8	16 - 1	134.2
17 - 1	28.4	17 - 1	134.5
20 - 1	27	20 - 1	138.5
22 - 1	29.8	22 - 1	132.8
23 - 1	28.5	23 - 1	140
24 - 1	28.5	24 - 1	135.9
25 - 1	28.3	25 - 1	138.2
26 - 1	28	26 - 1	133
27 - 1	29.8	27 - 1	138
29 - 1	29	29 - 1	136
30 - 1	28	30 - 1	138.2

Westgard	تفسیر	SGPT ( ALT )	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	28.5	SD	0.9	Exp. Date	1403/1/17	Lot. No	035

