

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.73	
Mean	28.11	Mean	134.53
SD	0.79	SD	2.33
CV	2.81	CV	1.73
Day-Run	Value	Day-Run	Value
1 - 1	27.2	1 - 1	135.2
2 - 1	27.5	2 - 1	132.7
4 - 1	28.1	4 - 1	136.5
5 - 1	28.4	5 - 1	136.4
6 - 1	27.3	6 - 1	132.8
7 - 1	27.7	7 - 1	139
8 - 1	29.3	8 - 1	132
11 - 1	28.6	11 - 1	133
12 - 1	27.5	12 - 1	136
13 - 1	28.6	13 - 1	130.7
14 - 1	28.9	14 - 1	135.7
18 - 1	29	18 - 1	136
19 - 1	30	19 - 1	137.5
20 - 1	26.8	20 - 1	132
21 - 1	28.4	21 - 1	130
22 - 1	27.6	22 - 1	135.9
26 - 1	27.6	26 - 1	135.9
27 - 1	27.9	27 - 1	134.9
28 - 1	27.9	28 - 1	133.8
29 - 1	27.8	29 - 1	134.6

Westgard	تفسیر	SGPT ( ALT )	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	27.9	SD	0.81	Exp. Date	1403/1/17	Lot. No	035

