

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 0.93	
Mean	39.21	Mean	144.19
SD	0.72	SD	1.35
CV	1.83	CV	0.93
Day-Run	Value	Day-Run	Value
3 - 1	40	3 - 1	145.8
4 - 1	39.5	4 - 1	142
5 - 1	37.9	5 - 1	141.5
6 - 1	39.4	6 - 1	143.2
8 - 1	39.1	8 - 1	142.7
9 - 1	39.8	9 - 1	143.7
10 - 1	39.7	10 - 1	145.7
12 - 1	39.4	12 - 1	146
13 - 1	40	13 - 1	145.4
15 - 1	40	15 - 1	143.1
16 - 1	39	16 - 1	145
17 - 1	38.8	17 - 1	145.2
20 - 1	39.2	20 - 1	144
22 - 1	39.4	22 - 1	145
22 - 2	38	22 - 2	142.8
23 - 1	39.6	23 - 1	145.6
24 - 1	39.4	24 - 1	144.7
25 - 1	37.6	25 - 1	145.6
26 - 1	39.6	26 - 1	143.6
27 - 1	38	27 - 1	144.8
29 - 1	39.8	29 - 1	144.2
30 - 1	39.4	30 - 1	142.6

Westgard	تفسیر	SGOT (AST)	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	38.73	SD	0.81	Exp. Date	1403/1/17	Lot. No	035

