

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 2.88	
Mean	4.22	Mean	5.03
SD	0.06	SD	0.14
CV	1.39	CV	2.88
Day-Run	Value	Day-Run	Value
1 - 1	4.2	1 - 1	5.1
2 - 1	4.2	2 - 1	5.2
4 - 1	4.3	4 - 1	4.8
5 - 1	4.2	5 - 1	5.2
6 - 1	4.2	6 - 1	5.1
7 - 1	4.3	7 - 1	5.2
8 - 1	4.3	8 - 1	4.9
11 - 1	4.1	11 - 1	4.9
12 - 1	4.2	12 - 1	5.2
13 - 1	4.2	13 - 1	4.8
14 - 1	4.2	14 - 1	4.9
18 - 1	4.2	18 - 1	5.1
19 - 1	4.2	19 - 1	5.2
20 - 1	4.3	20 - 1	4.9
21 - 1	4.2	21 - 1	5.1
22 - 1	4.2	22 - 1	5
26 - 1	4.2	26 - 1	5
27 - 1	4.2	27 - 1	5.1
28 - 1	4.3	28 - 1	4.8
29 - 1	4.1	29 - 1	5

Westgard	تفسیر	Albomin	آزمایش	BS-800 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	4.22	SD	0.1	Exp. Date	1403/1/17	Lot. No	035

