

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	80	
2 - 1	81	
4 - 1	78.9	
5 - 1	77.6	
6 - 1	83.1	
7 - 1	77.1	
8 - 1	77.8	
11 - 1	83.1	
12 - 1	80.7	
13 - 1	82	
14 - 1	79	
18 - 1	77.5	
19 - 1	80.8	
20 - 1	83	
21 - 1	78.6	
22 - 1	77.4	
26 - 1	79.4	
27 - 1	81	
28 - 1	79.2	
29 - 1	82.1	

Westgard	تفسیر	Amylase	آزمایش	BS-800 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	79.78	SD	2.2	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	79.97	SD	2.01	CV	2.52			

