

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	21.4	
3 - 1	21.45	
4 - 1	21.4	
5 - 1	21.1	
8 - 1	20.7	
9 - 1	21.31	
10 - 1	21.1	
12 - 1	20.98	
13 - 1	20.15	
15 - 1	21.2	
16 - 1	21.1	
17 - 1	20.9	
18 - 1	21.7	
19 - 1	21.3	
20 - 1	20.45	
22 - 1	21.4	
23 - 1	22	1:2S
24 - 1	21.3	
25 - 1	21.1	
26 - 1	20.7	
27 - 1	20.9	
29 - 1	19.9	1:2S

Westgard	تفسیر	BHCG	آزمایش	2000 ایمولایت XPI-	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	20.82	SD	0.45	Exp. Date	1403/2/10	Lot. No	36-02-028
Next Month	Mean	21.07	SD	0.48	CV	2.28			

