

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	21.45	1:2S
2 - 1	20.85	1:3S R4S
4 - 1	20.4	
5 - 1	20.9	
6 - 1	21.2	
7 - 1	21.1	
8 - 1	21.4	
9 - 1	20.85	
11 - 1	20.59	
12 - 1	20.62	
13 - 1	21	
14 - 1	20.5	
16 - 1	19.75	1:2S
18 - 1	20.12	
19 - 1	20.52	
20 - 1	21.2	
21 - 1	20.98	
22 - 1	21.25	

Westgard	تفسیر	BHCG	آزمایش	2000 ایمولایت XPI-	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	20.7	SD	0.37	Exp. Date	1402/4/25	Lot. No	36-02-028
Next Month	Mean	20.82	SD	0.45	CV	2.16			

