

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	0.99
SD	0.07
CV	7.55
Day-Run	Value
1 - 1	0.95
2 - 1	0.92
4 - 1	1.06
5 - 1	1.05
6 - 1	0.95
7 - 1	0.9
8 - 1	1.07
11 - 1	0.9
12 - 1	1.1
13 - 1	0.93
14 - 1	1.06
18 - 1	1.07
19 - 1	1.1
20 - 1	0.92
21 - 1	1.04
22 - 1	0.9
26 - 1	0.95
27 - 1	0.98
28 - 1	1.07
29 - 1	0.92

High - CV : 2.6	
Mean	4.27
SD	0.11
CV	2.6
Day-Run	Value
1 - 1	4.32
2 - 1	4.11
4 - 1	4.37
5 - 1	4.24
6 - 1	4.25
7 - 1	4.3
8 - 1	4.11
11 - 1	4.49
12 - 1	4.19
12 - 2	4.19
13 - 1	4.41
14 - 1	4.18
18 - 1	4.41
19 - 1	4.19
20 - 1	4.25
21 - 1	4.28
22 - 1	4.1
26 - 1	4.3
27 - 1	4.2
28 - 1	4.44
29 - 1	4.3

Westgard	تفسیر	Bili T	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	1	SD	0.06	Exp. Date	1403/1/17	Lot. No	035

