

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	91.7	
2 - 1	91.3	
4 - 1	91	1:2S
5 - 1	91.4	Cusum 3S
6 - 1	90.3	
7 - 1	91	
8 - 1	91	
9 - 1	91	
11 - 1	90.5	
12 - 1	90.8	
13 - 1	91.4	
14 - 1	91.6	
16 - 1	90.8	
18 - 1	90.2	
19 - 1	90.5	
20 - 1	90.5	
21 - 1	91.4	
22 - 1	91	
26 - 1	90.2	
27 - 1	91.1	
28 - 1	90.6	
29 - 1	90.9	
30 - 1	90.6	

Westgard	تفسیر	C3	آزمایش	BS-480 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	90.91	SD	0.44	Exp. Date	1402/9/14	Lot. No	641-01-030/2
Next Month	Mean	90.9	SD	0.43	CV	0.48			

