

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	18.1	
1 - 2	17.77	
4 - 1	18.1	
6 - 1	17.9	
7 - 1	17.9	
8 - 1	17.75	
9 - 1	17.9	
11 - 2	17.8	
12 - 1	17.9	
13 - 1	18	
14 - 1	18.11	
16 - 1	17.85	
18 - 1	17.9	
19 - 1	17.8	
20 - 1	17.9	
21 - 1	17.7	
22 - 1	17.8	
26 - 1	17.6	1:2S
27 - 1	17.9	
28 - 1	17.99	
29 - 1	17.8	
30 - 1	17.9	

Westgard	تفسیر	C4	آزمایش	BS-480 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	17.9	SD	0.13	Exp. Date	1402/9/14	Lot. No	641-01-030/2
Next Month	Mean	17.88	SD	0.13	CV	0.72			

