

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	50.48	
3 - 1	51.4	
4 - 1	51.1	
5 - 1	52.4	
6 - 1	50.14	
7 - 1	51.42	
8 - 1	50.48	
9 - 1	52.19	
10 - 1	51.85	
12 - 1	50.12	
13 - 1	49.7	
15 - 1	50.87	
16 - 1	52.6	
17 - 1	49.8	
18 - 1	51.4	
19 - 1	52.5	
20 - 1	47.4	
22 - 1	47.9	
23 - 1	49.1	
24 - 1	48.8	
25 - 1	51.7	
26 - 1	53.7	
27 - 1	47.6	
29 - 1	45.4	1:3S
30 - 1	47.1	1:2S

Westgard	تفسیر	CA 125	آزمایش	2000 ایمولایت XPI-	دستگاه	7	ماه	1402	سال
Control	Seronorm	Mean	50.5	SD	1.61	Exp. Date	1403/2/10	Lot. No	2011909
Next Month	Mean	50.49	SD	1.79	CV	3.55			

