

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	17.6	
2 - 1	16.5	
4 - 1	16.8	
5 - 1	17.1	Cusum 3S
7 - 1	17.4	
8 - 1	17.8	
9 - 1	17.6	
11 - 1	16.9	
12 - 1	18.12	
13 - 1	16.7	
14 - 1	16.5	
16 - 1	17.4	
18 - 1	17.8	
19 - 1	16.8	
20 - 1	16.93	
21 - 1	17.22	
22 - 1	18.2	
26 - 1	18.05	
27 - 1	17.68	
28 - 1	17.2	
29 - 1	16.5	
30 - 1	17.1	

Westgard	تفسیر	CEA	آزمایش	SIEMENS Centaur CP	دستگاه	6	ماه	1402	سال
Control	Seronorm	Mean	17.38	SD	0.48	Exp. Date	1403/6/19	Lot. No	365-02-028
Next Month	Mean	17.27	SD	0.54	CV	3.11			

