

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.39	
Mean	127.35	Mean	214.7
SD	2.16	SD	2.98
CV	1.7	CV	1.39
Day-Run	Value	Day-Run	Value
1 - 1	125	1 - 1	213
2 - 1	126	2 - 1	215
4 - 1	130	4 - 1	213
5 - 1	125	5 - 1	220
6 - 1	130	6 - 1	215
7 - 1	129	7 - 1	215
8 - 1	130	8 - 1	218
11 - 1	125	11 - 1	211
12 - 1	128	12 - 1	214
13 - 1	129	13 - 1	218
14 - 1	128	14 - 1	210
18 - 1	130	18 - 1	214
19 - 1	125	19 - 1	215
20 - 1	130	20 - 1	217
21 - 1	127	21 - 1	211
22 - 1	126	22 - 1	218
26 - 1	125	26 - 1	210
27 - 1	124	27 - 1	218
28 - 1	126	28 - 1	217
29 - 1	129	29 - 1	212

Westgard	تفسیر	CPK	آزمایش	BS-800 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	127.5	SD	2.06	Exp. Date	1403/1/17	Lot. No	035

