

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	47.5	
3 - 1	47.1	
4 - 1	47	
5 - 1	47.6	
6 - 1	47.8	
8 - 1	46.9	
9 - 1	47.8	
10 - 1	47.7	
12 - 1	47.4	
13 - 1	47.3	
15 - 1	46.9	
16 - 1	48	
17 - 1	47.2	
18 - 1	46.3	
19 - 1	46.6	
20 - 1	46.5	
22 - 1	47.6	
23 - 1	46.8	
24 - 1	47.2	
25 - 1	48	
26 - 1	46.6	
27 - 1	48	
29 - 1	48.3	

Westgard	تفسیر	CRP-quantitative	آزمایش	BS-480 اتوآنالایزر	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	47.29	SD	0.53	Exp. Date	1403/6/8	Lot. No	326-2-034/3
Next Month	Mean	47.31	SD	0.55	CV	1.17			

