

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	47.7	
2 - 1	47.3	
4 - 1	46.9	
5 - 1	48	
6 - 1	47.5	
7 - 1	47.5	
8 - 1	47.2	
9 - 1	47.5	
11 - 1	46.8	
12 - 1	47.5	
13 - 1	47.7	
14 - 1	48.3	
16 - 1	47.1	
18 - 1	47.3	
19 - 1	46.9	
20 - 1	47.3	
21 - 1	46.4	
22 - 1	46.3	
26 - 1	47.9	
27 - 1	47.3	
28 - 1	48	
29 - 1	46.6	
30 - 1	46.9	

Westgard	تفسیر	CRP-quantitative	آزمایش	BS-480 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	47.23	SD	0.55	Exp. Date	1403/6/8	Lot. No	326-2-034/3
Next Month	Mean	47.29	SD	0.53	CV	1.11			

