

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 2.44	
Mean	9.42	Mean	13.57
SD	0.14	SD	0.33
CV	1.52	CV	2.44
Day-Run	Value	Day-Run	Value
1 - 1	9.4	1 - 1	13.6
2 - 1	9.7	2 - 1	13.7
4 - 1	9.3	4 - 1	13.2
5 - 1	9.2	5 - 1	13.8
6 - 1	9.6	6 - 1	14.1
7 - 1	9.4	7 - 1	13.4
8 - 1	9.4	8 - 1	13.5
11 - 1	9.5	11 - 1	13.5
12 - 1	9.6	12 - 1	13.5
13 - 1	9.3	13 - 1	13
14 - 1	9.3	13 - 2	13
18 - 1	9.6	14 - 1	13.5
19 - 1	9.6	18 - 1	13.5
20 - 1	9.4	19 - 1	14
21 - 1	9.3	20 - 1	13.4
22 - 1	9.4	21 - 1	13.8
26 - 1	9.3	22 - 1	14.1
27 - 1	9.2	26 - 1	13.9
28 - 1	9.5	28 - 1	13.1
29 - 1	9.4	29 - 1	13.7

Westgard	تفسیر	Calcium	آزمایش	BS-800 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	9.41	SD	0.2	Exp. Date	1403/1/17	Lot. No	035

