

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	28.3	
3 - 1	30.6	
4 - 1	30.2	
5 - 1	29.5	
6 - 1	29.4	
8 - 1	28.3	
9 - 1	30.4	
10 - 1	30.1	
13 - 1	30.1	
15 - 1	29.6	
16 - 1	29	
17 - 1	28.8	
18 - 1	28.8	
19 - 1	29.8	
20 - 1	29.6	
22 - 1	29.8	
24 - 1	30.3	
25 - 1	29.7	
26 - 1	29.7	
27 - 1	30.2	
29 - 1	30	
30 - 1	29	

Westgard	تفسیر	Ceruloplasmin	آزمایش	BS-480 اتوآنالایزر	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	29.07	SD	0.92	Exp. Date	1406/9/14	Lot. No	641-01-030/2
Next Month	Mean	29.63	SD	0.64	CV	2.17			

