

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	30.2	
2 - 1	27.4	
4 - 1	30.3	
5 - 1	29.1	
6 - 1	29.8	
7 - 1	28.8	
8 - 1	29.2	
9 - 1	29.6	
11 - 1	27.8	
12 - 1	29.6	
13 - 1	29.8	
14 - 1	29.5	
16 - 1	27.7	1:2S
18 - 1	29.6	
19 - 1	27.2	
20 - 1	27.5	
21 - 1	29.7	
22 - 1	29.9	
26 - 1	29.2	
27 - 1	29.5	
28 - 1	29.4	
29 - 1	29	
30 - 1	28.9	

Westgard	تفسیر	Ceruloplasmin	آزمایش	BS-480 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	28.82	SD	0.87	Exp. Date	1406/9/14	Lot. No	641-01-030/2
Next Month	Mean	29.07	SD	0.92	CV	3.17			

