

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	175.48
SD	2.29
CV	1.31
Day-Run	Value
3 - 1	173
4 - 1	172
5 - 1	175
6 - 1	177
8 - 1	174
9 - 1	173
10 - 1	177
12 - 1	171
13 - 1	173
15 - 1	178
16 - 1	177
17 - 1	174
20 - 1	177
22 - 1	180
23 - 1	175
24 - 1	175
25 - 1	175
26 - 1	177
27 - 1	178
29 - 1	177
30 - 1	177

High - CV : 0.97	
Mean	263.19
SD	2.54
CV	0.97
Day-Run	Value
3 - 1	262
4 - 1	260
5 - 1	264
6 - 1	262
8 - 1	263
9 - 1	268
10 - 1	264
12 - 1	265
13 - 1	260
15 - 1	262
16 - 1	260
17 - 1	263
20 - 1	266
22 - 1	263
23 - 1	260
24 - 1	261
25 - 1	266
26 - 1	260
27 - 1	266
29 - 1	267
30 - 1	265

Westgard	تفسیر	Cholesterol	آزمایش	BS-800 اتوآنالایزر	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	175	SD	2.77	Exp. Date	1403/1/17	Lot. No	035

